



CATARACT EYE DISEASE DETECTION




MEMBER:

F2020376076 (Danish Jameel)



CATARACT:


CATARACT IS A COMMON EYE DISEASE THAT AFFECTS THE CLARITY OF THE LENS IN THE EYE, LEADING TO BLURRED VISION AND OTHER VISUAL IMPAIRMENTS. THE LENS OF THE EYE IS NORMALLY CLEAR AND HELPS TO FOCUS LIGHT ONTO THE RETINA AT THE BACK OF THE EYE, ALLOWING US TO SEE SHARP AND CLEAR IMAGES. HOWEVER, WITH A CATARACT, THE LENS BECOMES CLOUDED OR OPAQUE, OBSTRUCTING THE PASSAGE OF LIGHT AND CAUSING VISION PROBLEMS..



FACTOR BEHIND CATARACT:


CATARACTS CAN DEVELOP DUE TO A VARIETY OF FACTORS.

- AGE
- GENETICS
- TRAUMA
- DIABETES



THE SYMPTOMS OF CATARACTS CAN VARY, BUT COMMON SIGNS INCLUDE:

- **BLURRED OR HAZY VISION**
- **DIFFICULTY SEEING IN DIM OR LOW-LIGHT CONDITIONS**
- **INCREASED SENSITIVITY TO GLARE**
- **REDUCED COLOR PERCEPTION**
- **DOUBLE VISION IN ONE EYE**
- **FREQUENT CHANGES IN EYEGGLASS OR CONTACT LENS PRESCRIPTION**



About Datasets.


Ocular Disease Intelligent Recognition (ODIR) is a structured ophthalmic database of 5,000 patients with age, color fundus photographs from left and right eyes and doctors' diagnostic keywords from doctors.

This dataset is meant to represent “real-life” set of patient information collected by Shanggong Medical Technology Co., Ltd. from different hospitals/medical centers in China. In these institutions, fundus images are captured by various cameras in the market, such as Canon, Zeiss and Kowa, resulting into varied image resolutions.

Annotations were labeled by trained human readers with quality control management. They classify patient into eight labels including:


- Normal (N),
- Diabetes (D),
- Glaucoma (G),
- Cataract (C),
- Age related Macular Degeneration (A),
- Hypertension (H),
- Pathological Myopia (M),
- Other diseases/abnormalities (O)

DATA PREPROCESSING:

- **MOUNT GOOGLE DRIVE**
 - **READ A CSV**
 - **DEFINE A FUNCTION CALLED "HAS_CATARACT"**
 - **OBTAIN A SAMPLE OF NORMAL FUNDUS IMAGES FOR BOTH EYES TO BALANCE THE DATASET.**
 - **DEFINE THE DATASET DIRECTORY AND IMAGE SIZE (224)**
 - **SPLIT THE DATASET INTO FEATURES (X) AND LABELS (Y)**
- 
- A series of approximately 10-12 thin, black, parallel diagonal lines that originate from the bottom right corner of the slide and extend towards the top right, creating a sense of motion or a modern design element.

MODEL EXPLANATION:


1: RESNET-50:

- RESNET-50: RESNET-50 IS A DEEP CONVOLUTIONAL NEURAL NETWORK (CNN) ARCHITECTURE THAT CONSISTS OF 50 LAYERS.
 - THIS ALLOWS FOR EASIER TRAINING OF VERY DEEP NETWORKS AND HELPS TO ALLEVIATE THE DEGRADATION PROBLEM.
- 
- A series of approximately 10-12 parallel diagonal lines in the bottom right corner of the slide, extending from the bottom edge towards the right edge, creating a sense of movement or a modern design element.

MODEL EXPLANATION:

1: RESNET 101 -152:


- RESNET101 AND RESNET152: RESNET101 AND RESNET152 ARE VARIANTS OF THE RESNET ARCHITECTURE WITH 101 AND 152 LAYERS, RESPECTIVELY. THEY FOLLOW THE SAME RESIDUAL CONNECTION APPROACH AS RESNET-50 BUT HAVE A DEEPER ARCHITECTURE.
- THESE MODELS ARE EVEN MORE POWERFUL IN CAPTURING COMPLEX FEATURES



MODEL EXPLANATION:

3: VGG16 AND VGG19:

VGG16 AND VGG19 ARE CNN ARCHITECTURES DEVELOPED BY THE VISUAL GEOMETRY GROUP (VGG) AT THE UNIVERSITY OF OXFORD. THESE MODELS ARE CHARACTERIZED BY THEIR SIMPLICITY AND UNIFORMITY IN ARCHITECTURE. BOTH VGG16 AND VGG19 CONSIST OF MULTIPLE CONVOLUTIONAL LAYERS WITH SMALL 3X3 FILTERS FOLLOWED BY MAX POOLING LAYERS. THE MAIN DIFFERENCE BETWEEN VGG16 AND VGG19 IS THE NUMBER OF LAYERS, WITH VGG16 HAVING 16 LAYERS AND VGG19 HAVING 19 LAYERS.


A series of approximately 10-12 thin, black, parallel diagonal lines that originate from the bottom right corner and extend towards the top right corner of the slide, creating a sense of motion or a modern design element.

MODEL EXPLANATION:

4: DENSENET 169:

DENSENET-169 IS A CONVOLUTIONAL NEURAL NETWORK (CNN) ARCHITECTURE THAT BELONGS TO THE FAMILY OF DENSENET MODELS. IT WAS INTRODUCED BY HUANG ET AL. IN 2017 AS AN EXTENSION TO THE ORIGINAL DENSENET ARCHITECTURE.


THE MAIN BUILDING BLOCK OF DENSENET-169 IS THE "DENSE BLOCK," WHERE MULTIPLE CONVOLUTIONAL LAYERS ARE DENSELY CONNECTED TO FORM A BLOCK. WITHIN EACH DENSE BLOCK, THE OUTPUT FEATURE MAPS OF EACH LAYER ARE CONCATENATED TOGETHER AND SERVE AS INPUT TO SUBSEQUENT LAYERS WITHIN THE SAME BLOCK.

A series of approximately 10-12 parallel diagonal lines, slanted upwards from left to right, located in the bottom right corner of the image. They are black lines of uniform thickness.


MODEL EXPLANATION:

5: DENSENET 121:


DENSENET-121 IS ANOTHER CNN ARCHITECTURE THAT IS KNOWN FOR ITS DENSE CONNECTIVITY PATTERN. IN DENSENET, EACH LAYER IS DIRECTLY CONNECTED TO EVERY OTHER LAYER IN A FEED-FORWARD FASHION. THIS DENSE CONNECTIVITY HELPS IN IMPROVING FEATURE REUSE, REDUCING THE NUMBER OF PARAMETERS, AND ENHANCING GRADIENT FLOW. DENSENET-121 HAS 121 LAYERS AND HAS BEEN SHOWN TO ACHIEVE HIGH ACCURACY ON IMAGE CLASSIFICATION TASKS WHILE REQUIRING FEWER PARAMETERS COMPARED TO OTHER MODELS.

A series of parallel diagonal lines in the bottom right corner of the slide, extending from the bottom edge towards the right edge.


TECHNIQUES USED:

- EARLY STOPPING CALLBACK FUNCTION
 - SIGMOID FUNCTION
 - `MODEL.COMPILE(OPTIMIZER="ADAM",LOSS="BINARY_CROSS ENTROPY",METRICS=["ACCURACY"])`
 - FLATTEN IMAGE
 - DENSE LAYER
 - CROSS VALIDATION TECHNIQUE
- 
- A series of approximately 10 parallel diagonal lines extending from the bottom right corner towards the top right, creating a sense of motion or a modern design element.


RESULTS IN VISUALIZATIONS:




Precision Comparison for Each Model




Recall Comparison for Each Model




F1-score Comparison for Each Model




Accuracy Comparison for Each Model




Training Times of Different Models



Testing Accuracy of Different Models



Testing Loss of Different Models



CONCLUSION:



- The results provided include precision, recall, and F1-score values for the "Normal" and "Cataract" classes, as well as overall accuracy and additional metrics for the evaluation. The final conclusion can be summarized as follows:
- The model achieved an accuracy of 94.96% on the test dataset, indicating that it was able to correctly classify the majority of the images. The test loss, a measure of how well the model performed on the test set, was 1.4225.
- Considering the precision, recall, and F1-score for both the "Normal" and "Cataract" classes, the model achieved high scores for both classes, with values ranging from 0.92 to 0.98. This suggests that the model performed well in correctly identifying both normal and cataract cases.
- In conclusion, based on the provided results, the model appears to be effective in classifying normal and cataract cases with a high degree of accuracy.